

Centre for Human Metabolomics (CHM)

<b>Test:</b>	<b>Immunoreactive trypsinogen (IRT): BLOOD CARD SAMPLE [DBS]</b>				
<b>Test Code:</b>	4592 x1				
<b>Tariff (including VAT):</b>	R 453.27				
<b>Description:</b>	Assay, quantification and interpretation.?				
<b>Turnaround time:</b>	9 work days from receipt of sample at PLIEM laboratory				
<b>Transit stability / Sample viability:</b>	Keep in sealed paper envelope after dried according to requirements, send separate from other specimens and send within 2 days after collection. Viability: 1 month, kept in a dry place.				
<b>Comments:</b>	<ol style="list-style-type: none"> <li>1. Place dried blood card in sealed paper envelope and NOT in plastic bag.</li> <li>2. Blood card must NOT be placed in envelope before completely dry.</li> <li>3. Note: <b>Test can only be performed on patients younger than 24 days</b></li> </ol>				
<b>Sample required:</b>	1x Blood Card Sample – 4 complete circles				
<b>Method:</b>	Immuno assay				
<b>Reference ranges &amp; units:</b>	<p>Reference ranges in ng/ml:</p> <table border="1"> <thead> <tr> <th></th> <th>IRT</th> </tr> </thead> <tbody> <tr> <td>All ages</td> <td>&gt; 75</td> </tr> </tbody> </table>		IRT	All ages	> 75
	IRT				
All ages	> 75				
<b>Consultant/Scientist:</b>	Dr Marli Dercksen / Prof Chris Vorster				
<b>Telephone number:</b>	018 299 2302 / 018 299 4196				
<b>E-mail address:</b>	marli.dercksen@nwu.ac.za / chris.vorster@nwu.ac.za				
<b>Contact for results &amp; other enquiries:</b>	Ansie Mienie				
<b>Telephone number:</b>	018 299 2312 / 018 299 1815 / 018 285 2652 (leave message)				
<b>Fax number:</b>	018 299 2316				
<b>E-mail address:</b>	<a href="mailto:ansie.mienie@nwu.ac.za">ansie.mienie@nwu.ac.za</a>				
<b>Delivery address for samples:</b>	Att: Dr M Dercksen, School for Biochemistry, PLIEM LAB POTCH, Building F3, Laboratory G20 (ground floor), 11 Hoffmann street, North West University, Potchefstroom				